

# Wideband Handheld Confocal Raman Skin Analyzer **DermDive** Series



**2023 V1**


For customized projects please Contact us:

[info@simscop.com](mailto:info@simscop.com)

## **DermDive** - the first wideband handheld confocal Raman skin analyzer in the world

Confocal Raman Spectroscopy (CRS) is a groundbreaking optical technology used extensively for biological tissue analysis, prized for its noninvasive approach, pinpoint accuracy and comprehensive biochemical insights. Despite its advantages, traditional CRS systems are encumbered by limited detection, cumbersome designs, reduced flexibility, slow data collection, and escalated costs.

To counter these challenges, we have developed a compact, handheld CRS system, aptly named "**DermDive**". DermDive is the world first wideband handheld confocal Raman skin analyzer covering both finger-print region (FP, 450–1750  $\text{cm}^{-1}$ ) and high wavenumber region (HW, 2800–3800  $\text{cm}^{-1}$ ) with a handheld probe. **DermDive** has adjustable and user-defined measurement time suitable for various applications. It provides molecular information at different skin layers noninvasively, up to 200 $\mu\text{m}$ .




**DermDive** is designed to be user-friendly for skin diagnosis on various parts of human body, thanks to its flexibility.

## Features


- A harmonized integration of single or dual-wavelength lasers through an efficient (>90%) wavelength combiner, enabling near-instantaneous excitation.
- The groundbreaking Raman Spectra Separation Algorithm (PRSSA), a pioneer in spectrum separation, harnessing the maximum a posteriori probability (MAP) estimate. This offers an astounding >99% accuracy, as evidenced in both laboratory simulations and real-time human skin analyses.
- A notable reduction in data collection time, outpacing conventional CRS systems by over 50%.
  - ✓ Noninvasive
  - ✓ Fast
  - ✓ User friendly
  - ✓ Handheld
  - ✓ Depth profiling
  - ✓ Safe
  - ✓ Fingerprint region + high wavenumber region

## Light Path Diagram



## Available Raman Spectra Database

Natural moisturizing factors (NMF), water, ceramide, keratin, lactic acid, urea, urocanic acid at different depths




Spec.	DermDive
<b>Laser Source</b>	633nm 785nm
<b>Raman Range</b>	FP: 450 – 1750cm <sup>-1</sup> ; HW: 2800-3750cm <sup>-1</sup>
<b>Axial Resolution</b>	10um~20um
<b>Stratum Corneum Thickness*</b>	Up to 200um
<b>Sensitivity</b>	60dB
<b>Scanning time</b>	2-10s/depth
<b>Probe Head Size</b>	~12 cm × 6 cm × 5 cm
<b>Probe Head weight</b>	~800g
<b>System size</b>	~ Height: 1m, Width: 50cm, Length: 50cm
<b>System weight</b>	~ 40kg

### Testable Substances & Parameters

- Natural Moisturizing Factors(NMF):
  - Pyrrolidone carboxylic acid
  - Ornithine
  - Serine
  - Proline
  - Glycine
  - Histidine
  - Alanine
  - Ceramide
  - Keratin
  - Lactic acid
  - Urea
  - Urocanic acid
  - Water content
- 
- Stratum corneum thickness
  - Penetration depth of a certain chemical
  - Other chemicals of interest upon request


\*Stratum Corneum Thickness refers to the measurable skin thickness.

## Benchmarking with RiverD Benchtop CRS




Spectra of DermDive and RiverD at single depth and Various depths. Spectra qualities are comparable.

## Optical Performance Characterization




Axial resolutions and sensitivities with different optical components.

## Software



Software for automatic skin data acquisition at all depths



Software for on-the-spot skin constituents' analysis